# CO10BP Vitality™ Series Selectorized Leg Press / Calf Extension ∭



The reimagined Vitality<sup>™</sup> Series Selectorized Line provides a compact, value-priced conditioning circuit that fits the space and budget requirements for a wide range of facilities, from hospitality to multi-family residential to community centers. New contemporary styling allows you to enhance the Vitality Selectorized strength training experience with pieces from Vitality Benches & Racks and Resolute Strength<sup>™</sup> Selectorized.

The Vitality<sup>™</sup> Series Leg Press/Calf Extension features an articulating foot pad that stays vertical to mimic a squat movement and allows for full extension during calf exercise. Ratcheting seat back allows for desired start position. Exercisers can easily engage the add-on weight with a simple push of a lever to increase the work load.





GO.PRECOR.COM/C010BP © 2021 Precor Incorporated.

# CO10BP Vitality™ Series Selectorized Leg Press / Calf Extension ∭

### Walk-through Design

This convenient design allows exercisers to enter and exit on either side of the unit, a very beneficial feature when space is limited.

### Large Foot Platform

The large foot platform allows exercisers of all sizes to adjust foot placement as needed, and have the space to change positions for exercise variety.

### **Smooth Motion Path**

The design of the foot pad assembly ensures that there is a smooth and natural path of motion, much like that of a standing squat.

### **Helpful Guidance**

The conveniently located instructional placard provides step-bystep guidance on body position, movement and muscles worked.

## Specifications

Dimensions (L x W x H): 76 x 44 x 69 in/ 192 x 111 x 175 cm

Equipment Weight: 631 lb / 287 kg

Weight Stack: 295 lb / 134 kg

Add-on Weight: 5 lb / 2.5 kg

Frame and Finish: Heavy-duty steel tubing is fully welded in all structural areas. Powder coated frame.

Warranty: Visit www.precor.com for warranty terms.

# 

# Color Options\*

Frame



Black

phoistery

Black Pearl

\*Colors may vary.



