C015BP

Vitality™ Series Selectorized

Rear Delt / Pec Fly III



The reimagined Vitality™ Series Selectorized Line provides a compact, value-priced conditioning circuit that fits the space and budget requirements for a wide range of facilities, from hospitality to multi-family residential to community centers. New contemporary styling allows you to enhance the Vitality Selectorized strength training experience with pieces from Vitality Benches & Racks and Resolute Strength™ Selectorized.

The Vitality™ Series Rear Delt / Pec Fly features pivoting arms that adjust for exerciser arm length. The easy-to-adjust start position and dual hand positions enbale exercise variety. A long, narrow pad provides back support for Pec Fly and chest support for Rear Delt movements. Exercisers can easily engage the add-on weight with a simple push of a lever to increase the work load.

















CO15BP

Vitality™ Series Selectorized

Rear Delt / Pec Fly III

Adjustable Positions

An easy start position and dual hand positions provide variety for both the pec fly and rear delt movements.

Support for Both Movements

The long, narrow pad provides back support for Pec Fly exercise and chest support for Rear Delt movement.

Ensures Proper Positioning

To ensure proper positioning in both exercises with minimal adjustments, the arms automatically pivot when the exerciser enters, adjusting to various arm lengths and positions.

Helpful Guidance

The conveniently located instructional placard provides step-bystep guidance on body position, movement and muscles worked.

Specifications

Dimensions (L x W x H): 56 x 54 x 78 in / 141 x 136 x 199 cm

Equipment Weight: 431 lb / 196 kg Weight Stack: 240 lb / 109 kg Add-on Weight: 5 lb / 2.5 kg

Frame and Finish: Heavy-duty steel tubing is fully welded in all

structural areas. Powder coated frame.

Warranty: Visit www.precor.com for warranty terms.

Color Options*

Frame Black Pearl

Upholstery Black

*Colors may vary.



